

# Life Cycle Management

Imagination At Your Service



PATIENT  
SAFETY



IMPROVED  
EFFICIENCY



ENHANCED  
CLINICAL OUTCOMES



TAILORED  
SOLUTION

Every detail matters. Because every patient matters.



## AeroDR

## GE Healthcare and Konica Minolta partnering to enable you to move to digital

Digitalization of the hospital equipment is a major solution to improve work-flow and exams' quality while enhancing patient safety by decreasing the dose they receive.

But how can you achieve this transition without having to acquire an entire new system ?

With AeroDR, a digital panel solution that

- is adapted to your daily activity;
- brings you all the benefits of high technology;
- maximizes the TCO (Total Cost of Ownership) efficiency;
- is reliable and backed by GE Healthcare support.

GE Healthcare



Enjoy smooth daily operations with ergonomic, light and easy to handle panel, great wear & tear resistance and high performance battery.

Improve your workflow through it's universal connectivity.

Get an excellent Image quality with dramatic dose reduction through high technology panel.

Increase your patient throughput via robust design and full integration capabilities.



#### Key Specifications

Panel size	14x17" & 17x17" & 10x12"
Panel weight	2.9 kg (14x17") - 3.6 kg (17x17")
Scintillator technology	CsI (Cesium Iodid)
Full charge timing	30 minutes with the charging dock
Utilization time with one full charge	2 hours (120 images)
Battery life-time	7 years (10.000 full charging cycles)
Display of image after exposure	3 sec preview + 11 sec full
DQE (50% dose reduction vs CR)	65% @ 0 cycle/mm 51% @ 1 cycle/mm

For more information, please contact your local GE Healthcare Service representative.

#### About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter - great people and technologies taking on tough challenges. From medical imaging, software and IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

© 2014 General Electric Company - All rights reserved. General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information. GE and GE Monogram, InSite, Iling, OnWatch, AppsLinq are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company

JB18085AEb

