

Life Cycle Management

Imagination At Your Service



PATIENT
SAFETY



IMPROVED
EFFICIENCY



ENHANCED
CLINICAL OUTCOMES



TAILORED
SOLUTION

Every detail matters. Because every patient matters.



Mobile Radiographic Tables

The ideal partner for your C-arm

Flexibility and Speed are two key parameters for your C-arm. As your needs are constantly evolving, you also need a highly customizable solution.

Whether for angiography, hybrid OR procedures, or conventional endovascular:

Stille imagIQ2 is the ideal complement for your C-arm.



GE Healthcare

Fast and accurate positioning

through multidirectional and multi axial moves to avoid C-arm adjustments.

Allows to reduce dose without compromising on image quality

through good transparency rate.

Easy to adapt to your quickly evolving needs through a complete set of options for fully customized solution.

Features	Stille imagIQ2
patient Weight	225 (300kg w/load upg pack)
Table dimensions (L/W/H) in cm	234 / 77.2 / 71-109
Translucent area dimensions (L/W)	169.2 - 196 (headrest) / 46.3
Transparency (mm Al equivalent)	0.4
Movements in cm	True Free Float Long : 90/Lat: 25/Diag : 93.4 Vertical : 38 motorized
Tredelledburg	+/- 25°
Lateral (Iso) roll	+/- 15°
Power	AC 100-250V + DC battery

For more information, please contact your local GE Healthcare Service representative.

About GE Healthcare

GE Healthcare provides transformational medical technologies and services to meet the demand for increased access, enhanced quality and more affordable healthcare around the world. GE (NYSE: GE) works on things that matter - great people and technologies taking on tough challenges. From medical imaging, software and IT, patient monitoring and diagnostics to drug discovery, biopharmaceutical manufacturing technologies and performance improvement solutions, GE Healthcare helps medical professionals deliver great healthcare to their patients.

© 2014 General Electric Company - All rights reserved. General Electric Company reserves the right to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information. GE and GE Monogram, InSite, Iling, OnWatch, AppsLinq are trademarks of General Electric Company. GE Healthcare, a division of General Electric Company

JB18085AEed

